


Jenx Standz Size 1

Standz easily accommodates a range of abducted hip positioning for children between the ages of 1-9.

Supine/Prone/Upright/Abduction Standing Frame



Built in angle gauges for frame angle and individual leg abduction

Standz covers have been designed to wipe clean, they also contain an antimicrobial agent to aid infection control

Tool free adjustments - cleverly designed wing knobs offer a comfortable grip to enable effective clamping with minimal effort

Adjustable supports to accommodate knee flexion contractures and leg length discrepancies



Standard configuration

- Base with castors
- Hydraulic assisted tilt
- Pelvic and trunk support pads with strap
- Abducting legs
- Foot positioning plates
- Tray mount

The Standz is able to be used in the Prone or Supine configuration, depending on your child's needs



ACCESSORIES

Flat Pad Head Support
5240-4110-044



Multigrip Head Support
1 Link 5356-4601-000
2 Link 5356-4602-000
3 Link 5356-4603-000



Flat Shoulder Support
5240-2750-044



Multigrip Body Support
5240-2721-044



Pelvic Support Board
5240-1901-044



Back Pad
5240-2402-044



Standard Knee Support
5240-1750-044



Multi-Axis Knee Support
5240-1770-044



Sandals

Size 1 5356-6081-000
Size 2 5356-6082-000



Tray

Always recommended for Prone configuration
5240-7112-044

Tray Holder

5240-7112-044



Code	Angle Range	Age Range	Max User Height	Frame Width	Frame Length	Chest & Hip Pad Width	Knee Block Width	Height of Stander	Leg Abduction	Max User Weight	Foot Support Angle Dorsiflexion	Foot Support Angle Plantarflexion
5240-0005-044	0-90°	1-9 yrs	138cm	63.5cm	90cm	18-30cm	12cm	108cm	30°	45kg	10°	-25°